


HQ: Brumwell House, Westway Industrial Estate, Throckley, Newcastle Upon Tyne. NE15 9EW.

Protective Face Shield Visor – Reference 2002 Specification Sheet











Our Ref 2002 is a lightweight, clear face shield for personal use, that will create a barrier to help reduce droplets and sneezes.

USP

- The Protective Face Shield Visor Reference 2002 is lightweight in use and easy to wear
- The Face Shield can be cleaned with warm water and any kind of cleaning detergent
- Our visors are produced from Eco Friendly clear recycled APET.

Cleaning & Disinfection

- Our 2002 face shield can be cleaned with warm water and any kind of cleaning detergent

Application

- Personal use
- Help support essential business and NHS staff

Product Description

- A lightweight, clear face shield for personal use, that will create a barrier to help reduce droplets and sneezes. The parts used are removable, cleanable and replaceable.
- Scratched or damaged visors should be replaced

Features

- Very lightweight material
- Practical minimal design
- Reduction of droplets from coughs and sneezes

-

Advantages

- Extremely lightweight material
- Easily cleaned
- Reusable
- Long service life use

Storage Conditions

When stored in stacks, protect against heat, moisture and direct sunlight

Use & Handling

- Instructions on how to prepare and use the Face Shield have been provided
- Place your hand on the inside when putting on the Face Shield to avoid getting scratched by the edges

Period of obsolescence

Our 2002 is a single use product only to be used for no more than 12 hours

Cautions

- Do not apply much pressure when bending the Face Shield to avoid breaking it
- A warning that materials which may come into contact with the wearers skin could cause allergic reactions to susceptible individuals
- Eye-protectors against high speed particles worn over standard ophthalmic spectacles may transmit impacts, thus creating a hazard to the wearer

Packaging

- They are packed in a self-adhesive FDA food approved OPP bag. We pack our visors, straps, specification sheet and assembly instructions to customers' requirements.
- Shipped in standard shipping carton or jiffy bag depending upon quantity ordered

Product Material

- Protective Face Shield Visor Reference 2002 (456.1 x 240)
- 275 clear Micron Recycled APET
- Polyprop strap (700 Micron thick)

- *Disclaimer: Visor 2002 does not have a CE mark therefore customers purchase at their own risk UK Industrial Tapes Limited will not accept any liability.
- *Disclaimer: Whilst our Protective Face Shield Visor Reference 2002 has many applications it is not designed to be used with Covid-19 positive patients on its own without the use of other forms of PPE.

T: (0191) 269 7810 | F: (0191) 267 1028

E: sales@ukindustrialtapes.co.uk | sales@postapes.co.uk | sales@stix2.co.uk

Brumwell House, Westway Industrial Estate, Throckley, Newcastle Upon Tyne. NE15 9EW





