

## Protective Face Shield Visor – Reference 2002 Specification Sheet



**MADE IN  
BRITAIN** 



Our Ref 2002 is a lightweight, clear face shield for personal use, that will create a barrier to help reduce droplets and sneezes.

### **USP**

- The Protective Face Shield Visor – Reference 2002 is lightweight in use and easy to wear
- The Face Shield can be cleaned with warm water and any kind of cleaning detergent
- Our visors are produced from Eco Friendly clear recycled APET.

### **Cleaning & Disinfection**

- Our 2002 face shield can be cleaned with warm water and any kind of cleaning detergent

### **Application**

- Personal use
- Help support essential business and NHS staff

### **Product Description**

- A lightweight, clear face shield for personal use, that will create a barrier to help reduce droplets and sneezes. The parts used are removable, cleanable and replaceable.
- Scratched or damaged visors should be replaced

### **Features**

- Very lightweight material
- Practical minimal design
- Reduction of droplets from coughs and sneezes

### **Advantages**

- Extremely lightweight material
- Easily cleaned
- Reusable
- Long service life use

### **Storage Conditions**

- When stored in stacks, protect against heat, moisture and direct sunlight

### **Use & Handling**

- Instructions on how to prepare and use the Face Shield have been provided
- Place your hand on the inside when putting on the Face Shield to avoid getting scratched by the edges

### **Period of obsolescence**

- Our 2002 is a single use product only to be used for no more than 12 hours

### **Cautions**

- Do not apply much pressure when bending the Face Shield to avoid breaking it
- A warning that materials which may come into contact with the wearers skin could cause allergic reactions to susceptible individuals
- Eye-protectors against high speed particles worn over standard ophthalmic spectacles may transmit impacts, thus creating a hazard to the wearer

### **Packaging**

- They are packed in a self-adhesive FDA food approved OPP bag. We pack our visors, straps, specification sheet and assembly instructions to customers requirements.
- Shipped in standard shipping carton or jiffy bag depending upon quantity ordered

### **Product Material**

- Protective Face Shield Visor – Reference 2002 (456.1 x 240)
- 275 clear Micron Recycled APET
- Polyprop strap (700 Micron thick)

\*Disclaimer: Whilst our Protective Face Shield Visor - Reference 2002 has many applications it is not designed to be used with Covid -19 positive patients. However, we are currently working on a successor to our 2002 which is currently being tested in the clinical setting and CE approval has already been applied for.

T: (0191) 269 7810 | F: (0191) 267 1028

E: [sales@ukindustrialtapes.co.uk](mailto:sales@ukindustrialtapes.co.uk) | [sales@postapes.co.uk](mailto:sales@postapes.co.uk) | [sales@stix2.co.uk](mailto:sales@stix2.co.uk)

Brumwell House, Westway Industrial Estate, Throckley, Newcastle Upon Tyne, NE15 9EW

[www.ukindustrialtapes.co.uk](http://www.ukindustrialtapes.co.uk) | [www.postapes.co.uk](http://www.postapes.co.uk) | [www.stix2.co.uk](http://www.stix2.co.uk)



It is the buyer's responsibility to consider & assess the fitness for purpose of the goods, having regard to use to which the buyer intends to put such goods. The price is quoted for the quantity quoted. All prices quoted are for 30 days only. Standard credit terms are 30 days from the end of month unless agreed in writing otherwise. All orders are subject to this Company's terms & conditions of trading which can be seen on our website [ukindustrialtapes.co.uk](http://ukindustrialtapes.co.uk), [postapes.co.uk](http://postapes.co.uk) and [stix2.co.uk](http://stix2.co.uk)



POS Tapes & Stix2 are trading divisions to UK Industrial Tapes Ltd