


HQ: Brumwell House, Westway Industrial Estate, Throckley, Newcastle Upon Tyne. NE15 9EW.

C€ ₀₀₈₆

Protective Face Shield Visor – Reference 2005-Eco EN 166: 2002 Specification Sheet





Our Ref 2005-Eco is an eco-friendly biodegradable protective face visor. It is a lightweight, clear face shield for personal use, that will create a barrier to help reduce droplets and sneezes.

Here at UK Industrial Tapes we are fully committed to eco-friendly biodegradable* plastic products so we are very excited to be working with Breakdown Plastics Inc to bring you our new Eco-Friendly Range.

USP

- The Protective Face Shield Visor Reference 2005 is lightweight in use and easy to wear
- It has a Green Injection moulded Polyprop adjustable zip effect head strap
- Our Green Injection moulded head strap is very strong and durable and can be re-used every day
- The Face Shield can be cleaned with warm water and any kind of cleaning detergent
- Our visors are produced from Virgin PET with anti fog.

Cleaning, Disinfection and re - use

- Our 2005 face shield can be cleaned with warm water and any kind of cleaning detergent
- Our Green Injection Moulded head strap is very strong and durable and can be sanitized and cleaned every day ready for re use on a daily basis

Application

- Personal use
- Help support essential business and NHS staff

Product Description

- A lightweight, clear face shield for personal use, that will create a barrier to help reduce droplets and sneezes. The parts used are removable, cleanable and replaceable.
- Scratched or damaged visors should be replaced

Features

- Very lightweight material
- Practical minimal design
- Reduction of droplets from coughs and sneezes

Advantages

- Extremely lightweight material
- Easily cleaned
- Durable Re-usable green Injection Moulded head strap
- Long service life use

Breakdown Plastic Ince - The Science

Change Plastic for Good® offers a range of biodegradation-enhancing polymers and renewable resins that apply to a specific end-of-life scenario for plastic, rubber, foam, fabric and coated paper waste.

Breakdown Plastic Ince - In short

Ingredients in BDP™ accelerate the biodegradation of otherwise non-biodegradable plastics by providing a foothold for colonies of microorganisms. This foothold allows the microorganisms. This foothold allows the microorganisms to begin metabolizing the treated

BDP™ can be used in most plastics, synthetic fabric and yarn, rubber and foam

Storage Conditions

When stored in stacks, protect against heat, moisture and direct sunlight

Use & Handling

- Instructions on how to prepare and use the Face Shield have been provided
- Place your hand on the inside when putting on the Face Shield to avoid getting scratched by the edges

T: (0191) 269 7810 | F: (0191) 267 1028

E: sales@ukindustrialtapes.co.uk | sales@postapes.co.uk | sales@stix2.co.uk

Brumwell House, Westway Industrial Estate, Throckley, Newcastle Upon Tyne. NE15 9EW

www.ukindustrialtapes.co.uk | www.postapes.co.uk | www.stix2.co.uk











Period of obsolescence

- Our 2005 is a single use product only to be used for no more than 12 hours
- Replacement packs of 50 clear face shields can be purchased

Cautions

- Do not apply much pressure when bending the Face Shield to avoid breaking it
- A warning that materials which may come into contact with the wearers skin could cause allergic reactions to susceptible individuals
- Eye-protectors against high speed particles worn over standard ophthalmic spectacles may transmit impacts, thus creating a hazard to the wearer

Packaging

- They are packed in a self-adhesive FDA food approved OPP bag. We pack our visors, straps, specification sheet and assembly instructions to customers requirements.
- Shipped in standard shipping carton or jiffy bag depending upon quantity ordered

Replacement face Shields

We supply a pack or 50 replacement face shields which can be easily replaced onto our head strap after every shift

Product Material

- Protective Face Shield Visor Reference 2005 (210mm x 298mm)
- 275 clear Virgin PET with anti fog
- Green Injection moulded Polyprop adjustable zip tie effect head strap

*Disclaimer: Whilst our Protective Face Shield Visor - Reference 2005 has many applications it is not designed to be used with Covid-19 positive patients on its own without the use of other forms of PPE.

T: (0191) 269 7810 | F: (0191) 267 1028

E: sales@ukindustrialtapes.co.uk | sales@postapes.co.uk | sales@stix2.co.uk

Brumwell House, Westway Industrial Estate, Throckley, Newcastle Upon Tyne. NE15 9EW









OS TAPES